

Analysis conditions

- For liquid chromatography-mass spectrometry analysis, the user will provide a suitable chromatographic method for the separation of the sample components. Otherwise, the user should consult about method development and optimization conditions. The user will supply the columns and / or solvents and reagents if they are not included between the ones offered by the facility.
- Available columns at the facility:
 - 1) Zorbax SB-C18, 1.8 μm , 2.1 x 50 mm
 - 2) ACE C18, 3 μm , 150 x 2.1 mm
 - 3) Teknokroma C18 Mediterranea, 5 μm , 150 x 4.6 mm
 - 4) Teknokroma C18 Mediterranea, 5 μm , 250 x 4.6 mm
- Available solvent and reactants at the facility: water, acetonitrile, methanol, isopropanol, formic acid, acetic acid, triethylamine, NH_4OH , NH_4HCOO .